

PORTABLE OXYGEN CONCENTRATORS

AS OF LATEST FAA RULING

THIS DOCUMENT REVISED 2022

FAA APPROVED BATTERY DURATION	Setting 0.5	1.0	2.0	3.0	4.0	5.0	6.0
Caire AirSep Focus	--	--	1.5 hrs	--	--	--	--
Focus w/ Battery AirBelt	--	--	4.0 hrs	--	--	--	--
Caire AirSep Freestyle	--	6.5 hrs	3.5 hrs	3.0 hrs	--	--	--
FreeStyle Battery AirBelt	--	6.5 hrs	3.5 hrs	3.0 hrs	--	--	--
Caire AirSep FreeStyle 5	--	4.25 hrs	2.5 hrs	2.0 hrs	1.75 hrs	1.25 hrs	--
FreeStyle 5 Battery AirBelt	--	4.0 hrs	2.5 hrs	2.0 hrs	1.75 hrs	1.5 hrs	--
Belluscura XPLO2R (8 cell battery)	--	4hr 44min	4hr 20min	3hr 16min	2hr 20min	--	--
Caire FreeStyle Comfort (Single Battery)	--	8 hrs	4 hrs	3 hrs	2.25 hrs	2 hrs	--
Caire FreeStyle Comfort (Double Battery)	--	16 hrs	8 hrs	6 hrs	4.5 hrs	4 hrs	--
Inogen One G3 8 Cell (Also the OxyGo)	--	5.0 hrs	4.0hrs	3.0 hrs	2.0 hrs	1.5 hrs	--
G3 16 Cell Extended Battery	--	10.0 hrs	9.0 hrs	6.0 hrs	4.0 hrs	3.0 hrs	--
Inogen One G4 Single Battery	--	2hr 40min	2hr 15min	1hr 20min	--	--	--
G4 Double Battery (Also the OxyGo Fit)	--	5.0 hrs	4.5 hrs	2.5 hrs	--	--	--
Inogen One G5 Single Battery	--	6 hrs	4 hrs 45 min	3 hrs	2 hrs	1 hr 30 min	1 hr 15 min
G5 Double Battery (Also the OxyGo Next)	--	13 hrs	10 hrs	6 hrs 45 min	4 hrs 30 min	3 hrs 15 min	2 hrs 30 min
Inova Labs LifeChoice Activox	--	12.0 hrs	6.0 hrs	5.0 hrs	--	--	--
Activox External Battery +add	--	3.0 hrs	3.0 hrs	2.0 hrs	--	--	--
Inova Labs Activox Pro	--	12.5 hrs	7.25 hrs	6.25 hrs	--	--	--
Activox External Battery +add	--	3.0 hrs	3.0 hrs	3.0 hrs	--	--	--
Inova Labs Activox 4 L	--	10.25 hrs	8.25 hrs	5.0 hrs	4.0 hrs	--	--
Activox External Battery +add	--	3.0 hrs	3.0 hrs	3.0 hrs	3.0 hrs	--	--
Invacare Platinum Mobile	--	5.08 hrs	3.5 hrs	2.3 hrs	1.75 hrs	--	--
Invacare Platinum Mobile with 2nd Battery	--	10.17 hrs	7.0 hrs	4.67 hrs	3.5 hrs	--	--
Invacare XPO2	--	3.67 hrs	2.67 hrs	2.0 hrs	1.67 hrs	1.3 hrs	--
XPO2 Ext. Battery +add	--	3.5 hrs	2.4 hrs	2.0 hrs	1.5 hrs	1.0 hrs	--
O2 Concepts Oxlife Freedom 8 Cell	--	5 hrs	4 hrs	3 hrs	2 hrs	1.75 hrs	--
O2 Concepts Oxlife Freedom 12 Cell	--	10 hrs	8 hrs	6 hrs	4 hrs	2.5 hrs	--
Lifestyle P2 Rythm Healthcare	--	5 hr	4 hr	3 hr	2hr 14min	1hr 50min	--
Precision Medical Active5	--	6.5 hrs	4.3 hrs	2.7 hrs	2 hr	1.5 hrs	--
Precision Medical EasyPulse PM4150	--	4.7 hrs	3.4 hrs	2.7 hrs	2.0 hrs	1.6 hrs	--
EasyPulse Ext. Battery +add	--	4.7 hrs	3.4 hrs	2.7 hrs	2.0 hrs	1.6 hrs	--
Precision Medical EasyPulse PM4130	--	5.5 hrs	3.8 hrs	2.8 hrs	--	--	--
EasyPulse Ext. Battery +add	--	5.5 hrs	3.8 hrs	2.8 hrs	--	--	--
Philips Respironics SimplyGo Mini	--	6.0 hrs	4.5 hrs	3.0 hrs	2.25 hrs	2.0 hrs	--
SimplyGo Mini Double Battery	--	12.0 hrs	9.0 hrs	6.0 hrs	4.5 hrs	4.0 hrs	--
ResMed Mobi Internal	--	8.6 hrs	6.0 hrs	3.8 hrs	3.2 hrs	--	--
ResMed Mobi Ext Battery +add	--	4.5 hrs	3.0 hrs	1.6 hrs	1.6 hrs	--	--
DeVilbiss Healthcare iGo2	--	5 hrs	3.5 hrs	2.5 hrs	1.75 hrs	1.5 hrs	--
DeVilbiss Healthcare iGo	--	5.4 hrs	4.7 hrs	4.0 hrs	3.5 hrs	3.2 hrs	3.0 hrs
Continuous Flow	--	4.0 hrs	2.4 hrs	1.6 hrs	--	--	--
Precision Medical EasyPulse TOC	--	3.3 hrs	3.0 hrs	2.8 hrs	2.5 hrs	2.3 hrs	--
Continuous Flow 1/8th=5.0 hrs 1/4th=4.3 hrs	3.3 hrs	1.8 hrs	0.8 hrs	--	--	--	--
GCE Healthcare Zen-O	--	4.0 hrs	4.0 hrs	3.0 hrs	2.25 hrs	2.0 hrs	1.75 hrs
Zen-O w/ 2nd Battery	--	8.0 hrs	8.0 hrs	6.0 hrs	4.5 hrs	4.0 hrs	3.5 hrs
Continuous Flow Battery/with 2nd Battery	3.0/6.0hrs	1.75/3.5hrs	0.75/1.5hrs	--	--	--	--
Invacare Solo2	--	4.5 hrs	3.5 hrs	3.0 hrs	2.5 hrs	2.5 hrs	--
Continuous Flow	--	3.5 hrs	2.5 hrs	1.5 hrs	--	--	--
Oxlife Independence	3hr 8min	3.0 hrs	2hr 53min	2.0 hrs	1.5 hrs	1hr 8min	53 min
w/ 2nd Battery	6hr 15min	6.0 hrs	5hr 45min	4.0 hrs	3.0 hrs	2hr 15min	1hr 45min
Continuous Flow Battery/with 2nd Battery	2hr53min/5hr45min	2hr38min/5hr15min	1hr15min/2hr30min	45min/1hr30min	--	--	--
Philips Respironics SimplyGo	--	3.4 hrs	3.0 hrs	2.2 hrs	1.7 hrs	1.3 hrs	1.3 hrs
Continuous Flow	2.9 hrs	1.6 hrs	0.9 hrs	--	--	--	--
SeQual Eclipse	--	4.4 hrs	3.6 hrs	3.0 hrs	2.6 hrs	2.3 hrs	2.1 hrs
Continuous Flow	4.4 hrs	3.7 hrs	2.0 hrs	1.3 hrs	--	--	--

SeQual Eclipse 5

Pulse Delivery	Continuous Flow
Setting .5 NA	Setting .5 4.4 hrs
Setting 1.0 5.4 hrs	Setting 1.0 3.7 hrs
Setting 2.0 5.2 hrs	Setting 2.0 2.0 hrs
Setting 3.0 4.9 hrs	Setting 3.0 1.3 hrs
Setting 4.0 4.0 hrs	Pulse 128 mL 2.5 hrs
Setting 5.0 3.7 hrs	Pulse 160 mL 2.0 hrs
Setting 6.0 3.5 hrs	Pulse 192 mL 1.7 hrs



All times based on NEW batteries and are estimates as per each manufacturer. Each Airline has its own rules about flying with portable oxygen concentrators. Battery run times are shorter at higher settings.

1-800-346-3556

RENTALS & SALES 24HR SUPPORT INFO@AEROMEDIC.COM